

## Petri Dishes

### I. APPLICATIONS

1. Microbiology
2. Cell growth

### II. FEATURES

1. Materials: polystyrene
2. EO sterilized
3. Disposable
4. Glass-like clarity
5. Resistance temperature: 64°C~80°C

### III. PETRI DISH DESIGN

1. Available in 150 x 15mm, 90 x 15mm or 60 x 15mm
2. Triple-vented lid
3. Stackable

### IV. PACKAGING

1. Bulk, 500 pieces/case (90 x 15mm or 60 x 15mm)
2. Bulk, 200 pieces/case (150 x 15mm)

### V. DESCRIPTION

Item	Description	Packaging	Sterile
66-1515	150 x 15mm petri dish	Bulk, 200 pieces/case	Yes
66-1501	90 x 15mm petri dish	Bulk, 500 pieces/case	Yes
66-1560	60 x 15mm petri dish	Bulk, 500 pieces/case	Yes

### VI. STERILITY

1. Products are EO (Ethylene Oxide) sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.





### VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

**IX. PRODUCT PARAMETER**

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H
66-1501	Petri Dish	O.D. 89.3	O.D. 86.5	7.9	6.7	14.9	13.4	16.1	1.2
66-1560	Petri Dish	O.D. 57.4	O.D. 53.6	10.1	8.7	14.2	12.7	15.5	1.4
66-1515	Petri Dish	O.D. 150	O.D. 144	8.7	6.6	14.3	12.2	16.6	2.1



66-1501

66-1560

